

Cardiac Autonomic Neuropathy System Analyser - CANS 504



FEATURES:

- Analyses Sympathetic and Parasympathetic autonomic nervous function
- Simple array of (one minute) six tests completes the comprehensive study
- Fully automatic imported NIBP (Non-Invasive Blood Pressure) kit enhances result
- Display of raw ECG and Heart Rate Histogram supports the diagnosis
- Fully automated WINDOWS(PC operating system) based system
- Two full A4 page colour report with all patient details
- Resting ECG with resting HP
- Deep Breathing ECG, Expiration Inspiration ratio
- Supine to Standing ECG, R-R interval and 30:15 ratio
- Valsalva Maneuver ECG, R-R interval and Valsalva ratio
- Postural Hypertension resting BP, BP immediately after standing and 1 minute after standing, Fall in systolic BP assessed
- Sustained Handgrip BP before grip and BP after 5 minutes of sustained grip. Increase in diastolic BP assessed
- Final Impression consists of Parasympathetic function, Sympathetic function and Cardiac Autonomic Neuropathy function

DIABETIK FOOT CARE INDIA

(Madras Engineering Services)

Plot No. 38, 2/494, AGS Colony 3RD Main Road, Kottiwakkam, Chennai-600 041

Ph: 044 - 4356 4129, 3251 4129, 98400 93752

Email: mesmedi@vsnl.com : Website: www.diabeticfootcareindia.com